RΞCOVR LABS

RECOVR PRO

8-CHAMBER PRESSURE BOOTS

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ENGLISH - STARTS ON PAGE 2

Safety information	2
Package contents & product overview	3
LED touch control panel	4
Quick start guide	5
The six massage modes	6
Chamber exclusion & first session	7
After use, care & specifications	8
Simple troubleshooting	9

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ENGLISH

Recovr Pro - User Manual

IMPORTANT
Read this manual before first use. Keep it for future reference. This product is intended for wellness and recovery support and is not a substitute for medical diagnosis or treatment.

1. Safety information

- Do not turn the unit on before the boots are fitted and all air connections are secure.
- Do not unzip the boots while the unit is operating.
- Do not use the pump when it is wet. Keep the unit away from rain, water and excessive humidity.
- Do not use or store near hot surfaces, open flames or direct heat sources.
- Do not pull, crush, kink or sharply bend the hoses. Disconnect by holding the connector, not the hose.
- Do not open, repair or modify the pump. Contact an authorised service provider.
- Stop immediately if you experience pain, numbness, unusual swelling, breathing difficulty, dizziness or skin colour changes.

Consult a physician before use

Do not use unless cleared by a qualified healthcare professional if you have, or may have, acute inflammatory skin disease, arrhythmia, erysipelas, deep-vein thrombosis, pulmonary oedema, acute venous thrombosis, unstable hypertension, a pacemaker, impaired circulation or another condition affected by compression.

WARNING
The original manufacturer manual lists the conditions above but does not replace individual medical advice. If in doubt, do not use the product until a clinician confirms that compression therapy is appropriate.

Package contents

Check all parts before use

1	Main pump unit	1
2	1-to-2 air distributor	1
3	Air hoses	2
4	Power charger	1
5	8-chamber leg boots	2
6	Recovr Labs backpack / carry case	1

The picture is for reference only and not for making the original.

The supplied backpack is intended for storing and transporting the complete Recovr Pro set. Allow all parts to dry fully before packing.

Product overview

- Eight independently displayed compression chambers per boot.
- Six factory massage sequences (A-F).
- Adjustable treatment time: 10-90 minutes.
- Adjustable pressure: 30-260 mmHg in 5 mmHg increments (47 pressure levels).
- Rechargeable, portable pump with pause/resume and chamber exclusion controls.

LED touch control panel

Controls and indicators

The picture is for reference only and not for making the original.

TIME - / +	Reduce or increase session time.
PRESSURE - / +	Reduce or increase compression pressure.
POSITION	Move through chambers 1-8 to select a chamber.
CONFIRM	Exclude or restore the selected chamber.
MODE	Cycle through the six massage sequences.
ON/PAUSE	Start, pause or resume the session.
Battery icon	Shows the approximate remaining charge.
Chamber row 1-8	Shows active/selected chambers and inflation progress.

NOTE
The product photo shows the Recovr Labs control panel. Display segments may illuminate differently during operation.

Switching the pressure display

With the unit switched on, press and hold the CONFIRM button to switch between the current pressure during inflation and a steady display of the selected target pressure. Repeat to return to the previous display.

Quick start guide

Set up, select and start

Before using

Charge the pump before first use. While charging, the adapter indicator may show red; it changes to green when charging is complete. The battery icon on the pump indicates charge status.


- Place the pump on a stable, dry surface.
- Connect each air hose securely to the boots, then connect the main hose plug/distributor to the pump.
- Put on both boots comfortably. Smooth the material and close all zips/fasteners without trapping clothing or skin.
- Switch on the pump. The unit enters standby.
- Set time with TIME - / +. Range: 10-90 minutes; factory default: 30 minutes. Hold a key briefly for faster adjustment in 5-minute increments.
- Set pressure with PRESSURE - / +: 30-260 mmHg in 5 mmHg increments (47 levels). Factory default: 80 mmHg. Begin at a low, comfortable pressure.
- Select a massage sequence with MODE. The original manual specifies that mode changes are made while paused.
- Optional: press POSITION to select a chamber, then CONFIRM to exclude it from inflation. Repeat for additional chambers.
- Press ON/PAUSE to begin. Press again to pause or resume.

COMFORT CHECK
Compression should feel firm but comfortable. Stop the session if it becomes painful or if sensation or skin colour changes.

The six massage modes

How to read the inflation diagram

The boots contain eight air chambers, numbered from the foot toward the upper leg. In the diagram below, each small leg is the next step of the program. The black sections show which chambers are inflated at that moment. The sequence then repeats until the session ends or is paused.

A	Travelling wave: one chamber inflates at a time, moving from the foot toward the thigh. The previous chamber releases as the next one starts.
B	Build-up wave: compression starts at the foot and the lower chambers stay inflated while the next chambers are added, until the full leg is compressed.
C	Overlapping wave: neighbouring chambers work together and the compression zone moves upward in overlapping steps.
D	Alternating zones: separated groups of chambers inflate in a changing pattern, creating a more varied pulsing sensation.
E	Combination cycle: a travelling wave is followed by a build-up wave. This alternates a lighter pass with a fuller compression pass.
F	Extended combination: grouped and overlapping steps repeat from the foot upward, producing the most varied full-leg pattern shown in the guide.

CHOOSING A MODE
There is no single “best” mode. Start with A if you want the easiest pattern to understand, B for a fuller hold, or C/D for a more rhythmic feeling. E and F combine patterns. Comfort is more important than mode number or maximum pressure.

Excluding a chamber

- Pause the unit.
- Press POSITION until the required chamber number flashes.
- Press CONFIRM. The chamber light turns off, indicating that it will not inflate.
- Repeat if more chambers must be excluded. Verify the selection before restarting.

IMPORTANT
Chamber exclusion is a comfort feature, not medical programming. If a body area should not be compressed for clinical reasons, obtain professional advice before use.

A simple first session

- Sit or lie down comfortably and keep the pump uncovered on a firm surface.
- Leave all eight chambers active and select Mode A.
- Set 10-15 minutes and begin with low pressure.
- After the first full cycle, check that both legs feel comfortable and that the boots are not pinching.
- Increase pressure only a little at a time. Pause first if you want to change the mode.

What normal operation feels like

You will feel the pressure move or build according to the selected mode, followed by a release. A steady pump sound and changes in pressure are normal. Sharp pain, persistent numbness, burning, dizziness or breathing difficulty are not normal: stop immediately and remove the boots when they have deflated.

After use, care and specifications

Protect the hoses and keep the system clean

After use

- Pause and switch off the pump.
- Disconnect the pump-side connector by gripping the connector body. Never pull the hose.
- Deflate and remove the boots gently.
- Wipe down the boots and pump if required, allow them to dry, then store the set in the Recovr Labs backpack.

Cleaning and maintenance

- Switch off and unplug before cleaning.
- Do not wash or immerse the pump.
- Wipe with a dry or lightly damp soft cloth.
- Keep in a clean, dry place away from direct sunlight and heat.
- Inspect hoses, connectors, zips and cables regularly. Do not use damaged parts.

Pump specifications

Weight	1.5 kg
Dimensions	22.5 x 14.5 x 8.5 cm
Power consumption	49 W
Charger input	AC 100-240 V
Charger output	DC 12.6 V, 4 A
Battery	5200 mAh
Approx. charging time	2 hours
Approx. operating time	3-4 hours
Warranty stated in source manual	2 years

Specifications are reproduced from the supplied manufacturer manual. Actual labels supplied with the product and charger take precedence.

Simple troubleshooting

Checks you can safely make yourself

The pump does not switch on	Charge it fully. Check that the charger is firmly connected and that the wall socket works.
One or both boots do not inflate	Check every hose connection. Straighten any kinked hose. Make sure the session has started and the chambers have not been excluded.
The pressure feels too strong	Pause immediately and reduce the pressure. Choose a gentler mode. Stop completely if discomfort continues.
The pressure feels weak or uneven	Check for a loose connector, twisted hose or partly open zip. Confirm that both boots are fitted in a similar position.
The massage pumps in the wrong direction — from the thighs down to the feet	The hose connectors have been joined incorrectly. Disconnect them and reconnect each of the two hoses, making sure the small triangle on one side is precisely aligned with the matching mark on the opposite connector.
The pump makes a steady operating sound	A moderate pump sound is normal during inflation and release. Stop using it if the sound suddenly becomes very loud, irregular or is accompanied by heat or smell.
The unit stops by itself	The selected time may have ended or the battery may be low. Check the display and recharge if needed.

DO NOT OPEN THE PUMP
There are no user-serviceable parts inside. If the basic checks do not solve the problem, stop using the system and contact the seller or an authorised service provider.

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VETO S.A.